 NATIONAL
GEOGRAPHIC™

20X
MIKROSKOP
MICROSCOPE

- (DE)** Bedienungsanleitung
- (EN)** Operating Instructions
- (FR)** Mode d'emploi
- (NL)** Handleiding
- (IT)** Istruzioni per l'uso
- (ES)** Instrucciones de uso
- (PT)** Manual de utilização


Art. No. 9119100

MANUAL DOWNLOAD:

www.bresser.de/P9119100

EXPERIMENTS:


www.bresser.de/downloads

SERVICE AND WARRANTY:

www.bresser.de/warranty_terms

MICROSCOPE GUIDE:

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
MICROSCOPE FAQ:

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DE	Bedienungsanleitung	4
EN	Operating Instructions	6
FR	Mode d'emploi	8
NL	Handleiding	10
IT	Istruzioni per l'uso	12
ES	Instrucciones de uso	14
PT	Manual de utilização	16

**Garantie & Service / Warranty & Service / Garantie et Service /
Garantie et Service / Garantie & Service / Garanzia e assistenza /
Garantía y servicio / Garantia e Serviço** 18-19

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General Warnings

- **Choking hazard** – This product contains small parts that could be swallowed by children. This poses a choking hazard.
- Do not disassemble the device. In the event of a defect, please contact your dealer. The dealer will contact the Service Centre and can send the device in to be repaired, if necessary.

EC Declaration of Conformity

Bresser GmbH has issued a „Declaration of Conformity“ in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address:

www.bresser.de/download/9119100/CE/9119100_CE.pdf

Parts overview

1. WF 10x Eyepiece
2. 45° inclined monocular head
3. Microscope head
4. Upper locking screw (release for 360° rotation)
5. Objective lens
6. Locking Knob
7. Adjustable working distance up to 80 mm
8. Focus knob
9. Lower locking screw (release for 360° rotation)
10. Locking clamps
11. Microscope stage (with black and white object plate)
12. Transport case

1. What is a microscope?

A microscope contains two lens systems: the eyepiece and the objective. We're presenting these systems as one lens each so that the concept is easier to understand.

The lower lens (objective) produces a magnified image of the object. The picture, which you can't see, is magnified once more by the second lens (eyepiece, 1), which you can see as the 'microscope picture'.

2. Assembly and location

Before you start, choose an ideal location for using your microscope. It's important that you choose a spot with enough light for normal observation. Furthermore, it is recommended that you place the microscope on a stable surface, because a shaky surface will not lead to satisfactory results.

3. Normal observation

For normal observation, place the microscope in a bright location (near a window or desk lamp, for example).
Turn the focus knob (2) to the upper stop.

Now, place an object on the stage (4), directly under the objective (1). When you take a look through the eyepiece, you can see the magnified object. At this point, you still might see a slightly fuzzy picture. Adjust the image sharpness by slowly turning the focus knob (2).

You can rotate the microscope head 360° by loosening the upper locking screw. This makes it easier for you to work and observe the objects at the same time. To examine objects that do not fit between the work surface, loosen the lower locking screw and rotate the entire upper part of the microscope 360°.

4. Experiments

Use the following web link to find interesting experiments you can try out.

<http://www.bresser.de/downloads>

Notes on Cleaning

- Protect the device from dust and moisture.

Disposal

 Dispose of the packaging materials properly, according to their type, such as paper or cardboard. Contact your local waste-disposal service or environmental authority for information on the proper disposal.

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© 2018 National Geographic Partners LLC. All rights reserved.
NATIONAL GEOGRAPHIC and Yellow Border Design are trademarks of the National Geographic Society, used under license.