**DE** INSTRUCTION MANUAL  
**GB** BEDIENUNGSANLEITUNG  
**FR** MODE D'EMPLOI  
**NL** HANDLEIDING  
**IT** ISTRUZIONI PER L'OSU  
**ES** MANUAL DE INSTRUCCIONES  
**PT** MANUAL DE INSTRUÇÕES

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## General Information

### About this Instruction Manual

Please read the safety instructions in this manual carefully. To avoid damage to the unit and the risk of injury, only use this product as described in the manual. Keep the instruction manual handy so that you can easily look up information on all the functions.

#### **DANGER!**

You will find this symbol in front of every section of text which deals with the risk of severe injury or even death in the event of improper use.

#### **CAREFUL!**

You will find this symbol before every section of text that deals with the risk of minor to severe injuries resulting from improper use.

#### **NOTE!**

You will find this symbol in front of every section of text which deals with the risk of damage to property or the environment.

### Intended Use

This product is intended only for private use. It was developed for the magnified display of things in nature.

### General Warning

#### **DANGER!**

Tools with sharp edges and points are often used when working with this device. For this reason, store this device and all accessories and tools in a location that is out of the reach of children. There is a RISK OF INJURY!

#### **DANGER!**

This device contains electronic components which operate via a power source (power supply and/or batteries). Do not leave children unattended when using the device. Only use the device as described in the manual, otherwise you run the RISK of an ELECTRIC SHOCK.

#### **DANGER!**

Do not expose the device to high temperatures. Use only the supplied power adapter. Do not short circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short circuit, fires or even explosions.

#### **DANGER!**

Never bend, pinch, pull or run over the power and connecting cables or extensions/adapters. Protect the cables from sharp edges and heat. Before operation, check the device, cables and connections for damage. Never use a damaged unit or a unit with damaged power cables, etc. Damaged parts must be exchanged by an authorized service centre immediately.

#### **DANGER!**

Children should only use the device under adult supervision. Keep packaging material (plastic bags, rubber bands, etc.) out of the reach of children! There is a RISK OF CHOKING!

#### **CAREFUL!**

Children must not have access to the included chemicals and liquids. Do not drink the chemicals. Wash hands thoroughly with running water after use. In the event of contact with the eyes or mouth, rinse thoroughly with water. In the event of pain, contact a doctor immediately and take the substances with you.

### ! NOTE!

Do not disassemble the device. In the event of a defect, please contact your dealer. He will contact the Service Centre and can send the device in to be repaired, if necessary.

Do not expose the device to temperatures above 40°C.

### Fig. 1-2

- |                                   |                             |
|-----------------------------------|-----------------------------|
| ① Binocular sight                 | ⑩ Transmitted light unit    |
| ② Trinocular sight                | ⑪ Microscope arm clamp      |
| ③ Eyepiece                        | ⑫ Height adjustment clamp   |
| ④ Microscope table                | ⑬ Dioptre compensation      |
| ⑤ Power switch                    | ⑭ Barlow lens               |
| ⑥ Direct/transmitted light switch | ⑮ Bino-/trinocular toggling |
| ⑦ Focussing                       | ⑯ Adapter                   |
| ⑧ Zoom 1x - 4x                    | ⑰ Lens clamping screw       |
| ⑨ Direct light unit               | ⑱ Dimmer                    |

### PACKAGE CONTENTS

- Microscope with inbuilt lighting
- 2 „WF 10x“ eyepieces
- 2 „WF 20x“ eyepieces
- Blue colour filter
- Spare fuse
- Spare bulb
- Hood
- Black/white microscope plate
- Operating instructions

### HOW TO USE

1. First insert the power plug in a power socket
2. Then place a specimen you wish to examine on the microscope table and fix it in place with the clamps if necessary.
3. Press the On/Off switch (fig. 1, 5) on the back of the microscope base.
4. Turn the lighting setting knob to adjust the lighting as needed (fig. 1, 6). There are three specimen lighting options. Position „I“ activates only the direct light unit (fig. 2, 9). You can adjust this slightly vertically to give you the best possible lighting. Position „III“ activates only the transmitted light unit (fig. 2, 10). Position „II“ activates both lighting units. The OFF position switches all specimen lighting off. As your device has continuously controllable lighting (dimmer) (fig. 1, 18) optimal illumination of the object to be viewed is guaranteed.
5. You can then adjust the binocular sight to the gap between your eyes. Look through the binocular sight (fig. 1, 1) and adjust the space between the two eyepieces whilst doing so (fig. 1, 3) until you see only a single image.
6. The specimen will still not be wholly correctly focussed, though. Turn the focussing wheel (fig. 1, 7) until the image is really sharp.
7. Specimens of various sizes can be examined with the Advance ICD Microscope. Just change the lens unit height to adjust the microscope to the specific specimen. This is done by loosening the microscope arm clamp (fig. 2, 11) and adjusting the height to the specimen. When you're finished tighten the clamp again. To fix the lens unit in place undo the height limit clamp (fig. 2, 12) and push it until it is below the microscope arm. Then tighten the clamp again.
8. As eyesight differs dioptre compensation is built into your microscope. Adjust the rotatable eyepieces to their centre setting (fig. 2.13). Now look through the binocular sight. Close your right eye and use only the left to look through the eyepiece. Adjust until the specimen image is sharp. Now close the left eye and look through the eyepiece with the right eye. If the image is unclear turn the dioptre compensation (fig. 2, 13) until the image is in focus. The instrument is now adjusted for both your eyes.
9. To change magnification turn the zoom adjustment (fig. 1.8). This continuously adjusts magnification between 1x and 4x. Actual magnification is calculated by multiplying the zoom by the eyepiece value (10x or 20x).
10. To replace the lenses (fig. 1, 3) undo their fastening screws using a small screwdriver. Do not unscrew them too far. After replacing the lenses the screws can be carefully retightened.
11. To increase microscope magnification use the Barlow lens supplied (fig. 5). Screw the Barlow lens carefully into the mounting provided for it (fig 1, 14).

### Replacing the microscope table plate

In addition to the transparent plate a plate with a white and a black side is included. This plate is used to improve contrast.

**Using the trinocular sight.**

This is used for photography. If you use a Bresser MikrOkular it can easily be inserted in the trinocular extension. A mirror reflex camera can be added using the optional camera adapter and T2 ring.

Before you can use the trinocular sight the lever (fig. 1, 15) must be pulled out.

**CARE AND MAINTENANCE****Replacing the transmitted light unit bulb**

1. Use the ON/OFF switch (fig. 1, 5) to turn the microscope off and then remove the power plug.
2. Undo the clamping screw to open the housing and then remove the glass plate and blue filter. The bulb (fig. 4) and its mounting are now visible.
3. Carefully remove the old bulb from its socket.
4. Insert the new one carefully. Make sure the bulb is correctly inserted.
5. Clean it with a little alcohol to remove any fingerprints.
6. Then close the housing. The microscope is then again ready for use.

**Changing the bulb in the direct light unit**

1. Use the ON/OFF switch (fig. 1, 5) to switch the microscope off and unplug it.
2. To make the job easier remove the entire lens unit. This is done by undoing the clamping screw (fig. 1, 17). Then put the entire unit carefully aside.
3. Then carefully unscrew the filter mounting with filter. Make sure the filter doesn't fall out (fig. 4).
4. Carefully remove the old bulb from its socket.
5. Insert the new one carefully. Make sure the bulb is correctly inserted.
6. Screw the filter back on and fasten the entire unit back in place. The microscope is then again ready for use.

**Changing the fuse**

Two fuses are inbuilt to protect the microscope. Should one fail it can be changed as follows.

1. Turn the microscope off using the ON/OFF switch (fig. 1, 5) and then remove the plug.
2. Check which fuse needs replacing. To do so undo the fuse caps on the microscope back. This is best done using a coin.
3. Then check them both. The wire in either or both is often visibly broken. If so, then that fuse must be replaced. Alternatively, switch the fuses around to find out which is not functioning.
4. Your microscope is ready for use again once any defective fuse has been replaced.

**Keeping the lens clean**

Cover the microscope when not in use using the hood supplied to keep it clean.

**Your microscope is a precision instrument so please do not remove any part/s of it. Doing so may prevent it working properly.**

**TECHNICAL DATA**

|                           |                          |
|---------------------------|--------------------------|
| Working gap               | max: 200 mm / min: 25 mm |
| Maximum specimen height   | approximately 120 mm     |
| Eye separation adjustment | 52-75 mm                 |
| Lighting                  | 12V 10W halogen bulb     |
| Lighting (direct)         | 12V 10W halogen bulb     |
| Mains power needed        | 230V 50HZ                |

Magnification and field of view:

| Eyepiece:                       | Magnification: | Field of view:       |
|---------------------------------|----------------|----------------------|
| 10x wide-angle                  | 10 - 40x       | approx. 22 - 5,5 mm  |
| 10x wide-angle with Barlow lens | 20 - 80x       | approx. 11 - 2,5 mm  |
| 20x wide-angle                  | 20 - 80x       | approx. 11 - 2,5 mm  |
| 20x wide-angle with Barlow lens | 40 - 160x      | approx. 5,5 - 1,1 mm |

**DISPOSAL/RECYCLING**

Dispose of the packaging materials properly, according to their type, such as paper or cardboard. Contact your local waste-disposal service or environmental authority for information on the proper disposal.

Do not dispose of electronic devices in the household garbage!

As per Directive 2012/19/EU of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner.

**EC DECLARATION OF CONFORMITY**

Bresser GmbH has issued a „Declaration of Conformity“ in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address:

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**UKCA DECLARATION OF CONFORMITY**

Bresser GmbH has issued a „Declaration of Conformity“ in accordance with applicable guidelines and corresponding standards. The full text of the UKCA declaration of conformity is available at the following internet address:

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Great Britain

**WARRANTY**

The regular warranty period is 2 years and begins on the day of purchase. To benefit from an extended voluntary warranty period as stated on the gift box, registration on our website is required.

You can consult the full warranty terms as well as information on extending the warranty period and details of our services at [www.bresser.de/warranty\\_terms](http://www.bresser.de/warranty_terms).

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